

[Home](#) >> [Junior & Childrens Range](#) >> [Junior and childrens walkers](#) >> Nimbo Lightweight Posterior Posture Walker

Product Information for [Nimbo Lightweight Posterior Posture Walker](#) on www.discountmobilityshop.com

Product Details

Title: Junior walking frame, Childrens walkers, Pediatric Walking frame, disabled child walker, kids walking frame

Product Name: [Nimbo Lightweight Posterior Posture Walker](#)

Category: [Junior and childrens walkers](#)

Model: KA1200N KA2200N KA3200N KA4200N

Description:

Nimbo Lightweight Posterior Posture Walker from Drive **RRP £165** Supports the user from behind, promoting upright posture and trunk extension. Lightweight aluminium, height adjustable complemented with one directional soft rubber rear wheel stopping the walker from rolling backwards. Front wheels swivel to enhance manoeuvrability and are skid resistant, ideal for indoor and outdoor use. With latex-free handgrips with flanged ends. Many therapists prefer posterior (reverse) walkers for children with cerebral palsy and other neurological conditions. Available in four sizes. KA1200N TYKE, GOLDROD YELLOW KA2200 JUNIOR, LIME GREEN KA3200N YOUTH, CORNFLOWER BLUE KA4200N YOUNG ADULT, MIDNIGHT BLUE The Nimbo Pelvic Stabilizer The pelvic stabilizer is adjustable in width/height and features back & lateral pads and a support belt. It is recommended for users who have an asymmetrical posture when standing or walking and therefore, cannot stay centered in the walker. KA1050 for use with all Nimbo Posture Walkers The Nimbo Seat The Nimbo seat mounts on to the standard walker. It folds up for standing and walking, plus folds down for convenience. KA1285 - For use with KA1200 KA2285 - For use with KA2200 KA3285 - For use with KA3200 KA4285 - For use with KA4200 Nimbo Forearm Platform Constructed of lightweight aluminium, the Nimbo forearm platforms are height and depth adjustable. The handgrips extend forward or back to the position needed, supinate, pronate and are angle adjustable. Each side is individually adjustable and allows custom positioning for each user. The contoured armrest is molded in one piece with a flexible edge for added - comfort and has Velcro straps to prevent arms from slipping. They can be mounted to the mounting bracket either towards the inside or outside of the walker for lateral adjustment. Spacers are included for additional lateral range. KA1035 - Forearm Platform, for use with all Nimbo Posture Walkers KA1035S - Mounting Bracket, for use with all Nimbo Posture Walkers Order on line or call Us CALL FREEPHONE : 08000 430357

Weight: 6.00 kg

Bonus Points: 93 points

Product Variables

Which Walker ? : Tyke Yellow KA1200N, Junior Lime Green KA2200N, Youth Blue KA3200N, Young Adult KA4200N

Extras : Forearm Platform, Seat 1285 for KA1200, Seat 2285 for KA2200, Seat 3285 for KA3200, Seat 4285 for KA4200, Pelvic Stabilizer